HEINE Oftalmoscopios Directos

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Electromagnetic disturbances – requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
	The supply voltage quality should be that of a typical hospital environment.
	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Performance features of the ME system that have been determined to be essential to the performance	None
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	
	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
Warning	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the considered device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

A list of all cables, transducers and other accessories that are relevant for the EMC compliance

Test	Compliance
RF emissions CISPR11	Group 1 Class B
Conducted emissions (EN 55011 CISPR 11)*	
Radiated emissions (EN 55011 CISPR 11)	
Harmonic current emissions (IEC 61000-3-2)*	Passed
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*	
Immunity	See attached immunity test levels

Immunity test levels	Test level	
Test	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact discharge: ±8 kV Air discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m 80 – 2700 MHz 80 % AM at 1kHz	
Electrical fast transients bursts (IEC 61000-4-4)*	± 2 kV 100 kHz repetition frequency	Direct ophthalmoscope with a 2,5 V bulb: not applicable, because of the internally powered mode. Direct ophthalmoscope with a 3,5 V bulb: compliance test levels are identically to the IEC 60601-1-2 test levels.
Surges (IEC 61000-4-5)*	± 0.5 kV, ± 1 kV	Direct ophthalmoscope with a 2,5 V bulb: not applicable, because of the internally powered mode. Direct ophthalmoscope with a 3,5 V bulb: compliance test levels are identically to the IEC 60601-1-2 test levels.
Conducted disturbances induced by RF fields (IEC 61000-4-6)*	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Direct ophthalmoscope with a 2,5 V bulb: not applicable, because of the internally powered mode. Direct ophthalmoscope with a 3,5 V bulb: compliance test levels are identically to the IEC 60601-1-2 test levels.
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)		385 MHz; Pulse modulation: 18 Hz; 27 V/m 450 MHz, FM: ±5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse modulation: 217 Hz; 28 V/m 2450 MHz; Pulse modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse modulation: 217 Hz; 9 V/m
Power frequency magnetic fields (IEC 61000-4-8)		30 A/m; 50Hz or 60 Hz
Voltage dips (IEC 61000-4-11)*	0 % U _T ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _T ; 1 cycle and 70 % U _T ; 25 30 cycles single phase: at 0°	Direct ophthalmoscope with a 2,5 V bulb: not applicable, because of the internally powered mode. Direct ophthalmoscope with a 3,5 V bulb: compliance test levels are identically to the IEC 60601-1-2 test levels.
Short interruptions (IEC 61000-4-11)*	0 % U _T ; 250 300 cycle	

*n | a : „Not applicable“ in the internally powered mode

Technical specification	
Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Nominal voltage	3.0 V – 3.7 V
Nominal current	3

Instruction for use HEINE direct ophthalmoscopes

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